



Recombinant Schizosaccharomyces pombe DNA primase small subunit (pri1)

Product Code	CSB-YP018680SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O14215
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	<p>MTVQIDELDD KDLDEIIANG TLDGAKQGAV DSETMIQYYR HLFPWKYLFQ WLNHGPPVVTN DFAHREFAFT LPNDAYIRYL SFSNWEELKK EALNLCPSRF EVGPVYSANP RDRKTIRKST FHPLKKELVF DIDMTDYDDV RTCCSKTNIC EKCWPFITIA VQVLDICFHE DFGFKHILWV YSGRRGIHAW ICDEIACSLD DRSRMIASY LQVVVGPNPQG GVRLINNLKR PLPHLTRSL NILKSAFVKI VLEDQDPWAS KEGAENLLKL LPDKDLASAL RKKWEVDPER SSKNKWSDID TVLASGSIAS ISPSVIAIAK QDIVLTYLYP RLDVEVSRHL NHLLKSPFCV HPGTSRVCVP IDIERMDSFN PLKVPTVNDL LQELDKNSQN DNGHGPTMET NTTENQKDNA RGQSNKGHGF STSLNPYTLY FKSFSSQLFK ETVGNKRKHE NLEF</p>
Source	Yeast
Target Names	pri1
Protein Names	Recommended name: DNA primase small subunit EC= 2.7.7.-
Expression Region	1-454
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.