



MFGE8 Antibody, Biotin conjugated

Product Code	CSB-PA07974D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08431
Immunogen	Recombinant Human Lactadherin protein (25-334AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lactadherin (Breast epithelial antigen BA46) (HMFG) (MFGM) (Milk fat globule-EGF factor 8) (MFG-E8) (SED1) [Cleaved into: Lactadherin short form; Medin], MFGE8
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	MFGE8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.