




Bdnf Antibody

Product Code	CSB-PA16639A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23363
Immunogen	Recombinant Rat Brain-derived neurotrophic factor protein (137-237AA)
Raised In	Rabbit
Species Reactivity	Rat, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Brain-derived neurotrophic factor (BDNF), Bdnf
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Neuroscience
Target Names	Bdnf

Image



Western blot
All lanes: Bdnf antibody at 2µg/ml + U251 whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 29 kDa
Observed band size: 29 kDa