




Timp1 Monoclonal Antibody

Product Code	CSB-MA167191A0m
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30120
Immunogen	Recombinant Rat Timp1
Raised In	mouse
Species Reactivity	Rat
Specificity	specific for Rat Timp1 denatured and native forms
Tested Applications	ELISA,IHC;Recommended dilution:IHC:1:50-1:500
Relevance	Complexes with metalloproteinases (such as collagenases) and irreversibly inactivates them by binding to their catalytic zinc cofactor. Also stimulates steroidogenesis by Leydig and ovarian granuloma cells; procathepsin L is required for maximal activity.
Form	liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%,Protein G purified
Isotype	IgG1
Clonality	monoclonal
Alias	Timp-1
Product Type	Monoclonal Antibody
Immunogen Species	Rattus norvegicus (Rat)
Target Names	Timp1
Clone No.	4A2B1

Image



Immunohistochemistry of paraffin-embedded Rat prostate cancer using CSB-MA167191A0m at dilution of 1:100.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.