



Ccl25 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA09189A0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O35903 |
| Immunogen | Recombinant Mouse C-C motif chemokine 25 protein (25-144AA) |
| Raised In | Rabbit |
| Species Reactivity | Mouse |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | C-C motif chemokine 25 (Chemokine TECK) (Small-inducible cytokine A25) (Thymus-expressed chemokine), Ccl25, Scya25 Teck |
| Immunogen Species | Mus musculus (Mouse) |
| Research Area | Others |
| Target Names | Ccl25 |