



Recombinant Human Actin-related protein 2/3 complex subunit 4 (ARPC4)

Product Code	CSB-MP002130HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P59998
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	TATLRPYLS AVRATLQAAL CLENFSSQVV ERHNKPEVEV RSSKELLQPV VTISRNEKEK VLIEGSINSV RVSIAVKQAD EIEKILCHKF MRFMMMRAEN FFILRRKPVE GYDISFLITN FHTEQMYKHK LVDFVIHFME EIDKEISEMK LSVNARARIV AEEFLKNF
Source	Mammalian cell
Target Names	ARPC4
Protein Names	Recommended name: Actin-related protein 2/3 complex subunit 4 Alternative name(s): Arp2/3 complex 20 kDa subunit Short name= p20-ARC
Expression Region	2-168
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes one of seven subunits of the human Arp2/3 protein complex. The Arp2/3 protein complex has been implicated in the control of actin polymerization in cells and has been conserved through evolution. The exact role of This protein, the p20 subunit, has yet to be determined. Three transcript variants encoding two distinct isoforms have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.