



Cxcl3 Antibody

Product Code	CSB-PA09239A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6W5C0
Immunogen	Recombinant Mouse C-X-C motif chemokine 3 protein (32-100AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-X-C motif chemokine 3 (Dendritic cell inflammatory protein 1), Cxcl3, Dcip1 Gm1960
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Cxcl3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.