



Cxcl14 Antibody

Product Code	CSB-PA09319A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9WUQ5
Immunogen	Recombinant Mouse C-X-C motif chemokine 14 protein (23-99AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-X-C motif chemokine 14 (B-cell and monocyte-activating chemokine) (Chemokine BRAK) (Kidney-expressed chemokine CXC) (MIP-2G) (Small-inducible cytokine B14), Cxcl14, Bmac Kec Ks1 Mip2g Scyb14
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Cxcl14