



IdcA Antibody

Product Code	CSB-PA303469HA01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P76008
Immunogen	Recombinant Escherichia coli Murein tetrapeptide carboxypeptidase protein (1-304AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.1mol/L NaCl, 15mmol/L NaN3, PH7.2
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Murein tetrapeptide carboxypeptidase (EC 3.4.17.13) (LD-carboxypeptidase A) (Muramoyltetrapeptide carboxypeptidase), IdcA, ycgQ
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	IdcA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.