



OIH Antibody, FITC conjugated

Product Code	CSB-PA318020LC01CH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10184
Immunogen	Recombinant Gallus gallus Ovoinhibitor protein (26-472AA)
Raised In	Rabbit
Species Reactivity	Gallus gallus
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ovoinhibitor, OIH
Immunogen Species	Gallus gallus (Chicken)
Research Area	Others
Target Names	OIH