



OIH Antibody

Product Code	CSB-PA318020LA01CH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10184
Immunogen	Recombinant Gallus gallus Ovoinhibitor protein (26-472AA)
Raised In	Rabbit
Species Reactivity	Gallus gallus
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ovoinhibitor, OIH
Immunogen Species	Gallus gallus (Chicken)
Research Area	Others
Target Names	OIH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.