



Efna5 Antibody, FITC conjugated

Product Code	CSB-PA007464LC01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P97605
Immunogen	Recombinant Rat Ephrin-A5 protein (21-203AA)
Raised In	Rabbit
Species Reactivity	Rat
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ephrin-A5 (AL-1) (EPH-related receptor tyrosine kinase ligand 7) (LERK-7), Efna5, Eplg7 Lerk7
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Neuroscience
Target Names	Efna5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.