



Efna5 Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA007464LC01RA |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P97605 |
| Immunogen | Recombinant Rat Ephrin-A5 protein (21-203AA) |
| Raised In | Rabbit |
| Species Reactivity | Rat |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Ephrin-A5 (AL-1) (EPH-related receptor tyrosine kinase ligand 7) (LERK-7), Efna5, Eplg7 Lerk7 |
| Immunogen Species | Rattus norvegicus (Rat) |
| Research Area | Neuroscience |
| Target Names | Efna5 |