



HLA-DRA Antibody, FITC conjugated

Product Code	CSB-PA17939C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01903
Immunogen	Recombinant Human HLA class II histocompatibility antigen, DR alpha chain protein (28-254AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HLA class II histocompatibility antigen, DR alpha chain (MHC class II antigen DRA), HLA-DRA, HLA-DRA1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	HLA-DRA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.