




HLA-DRA Antibody

Product Code	CSB-PA17939A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01903
Immunogen	Recombinant Human HLA class II histocompatibility antigen, DR alpha chain protein (28-254AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HLA class II histocompatibility antigen, DR alpha chain (MHC class II antigen DRA), HLA-DRA, HLA-DRA1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	HLA-DRA

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA17939A0Rb at dilution of 1:20

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.