



Recombinant Eliomys quercinus Phosphatidylcholine-sterol acyltransferase (LCAT)

Product Code	CSB-EP012783EKD-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O35573
Product Type	Recombinant Protein
Immunogen Species	Eliomys quercinus (Garden dormouse)
Purity	>85% (SDS-PAGE)
Sequence	EDCFTIWLDL NIFLSLGVDC WIDNTRVIYN RSSGYMSNAP GVQIRVPGFG KTYSIEYLDD NKLAGYMHTL VQNLVNNAYV RDETVRAPPY DWRLEPRHQE EYYLKLGLV EEMYATYGKP VFLIGHSLGF CHLLYFLLLQ PQGIPIMSSI KLVEEQRITT TSPWMFPSHQ VWPEDHVFIS TPNFNFTFSD FQRFFADLHF EDGWYMWLQS RDLLAGLPAP GVEVYCLYGV GLPTPHTYMY DHGFPYDTPV GIIYEDGDDT VTTHSIELCS HWQGRQPQPV HLLPLRGTQH LNMVFSNKT
Source	E.coli
Target Names	LCAT
Protein Names	Recommended name: Phosphatidylcholine-sterol acyltransferase EC= 2.3.1.43 Alternative name(s): Lecithin-cholesterol acyltransferase Phospholipid- cholesterol acyltransferase
Expression Region	1-299
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.