



# Recombinant Tatera kempfi gambiana Phosphatidylcholine-sterol acyltransferase (LCAT)

<b>Product Code</b>	CSB-EP012783TLK-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O35840
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Gerbilliscus gambianus (Gambian gerbil) (Gerbilliscus kempfi gambiana)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	TIWLDLNMFL SLGVDCWIDN TRVVYNRSSG RVSNAPGVEI RVPFGFKTYS VEYLDDNKLA EYMHTLVQNL VNNGYVRDET VRAAPYDWRL EPSQQDDYYQ KLAGLIEEMY AAYGKPVFLI GHSLGCLHVL YFLLRQGIPI MSSIKLREEQ RITTTSPWMF PDRDVPEDH VFISTPEFNY TGQDFERFFS DLHFEEGWYM WLQSRDLLAG LPAPGVDVYC LYGVGLPTPH TYIYDHNFPY KDPVAALYED GDDTVATRST ELCGQWQGRQ SQPVHLLPMN GTEHLNMVFS NKT
<b>Source</b>	E.coli
<b>Target Names</b>	LCAT
<b>Protein Names</b>	Recommended name: Phosphatidylcholine-sterol acyltransferase EC= 2.3.1.43 Alternative name(s): Lecithin-cholesterol acyltransferase Phospholipid- cholesterol acyltransferase
<b>Expression Region</b>	1-293
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.