



Recombinant Tatera kempfi gambiana Phosphatidylcholine-sterol acyltransferase (LCAT)

Product Code	CSB-YP012783TLK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O35840
Product Type	Recombinant Protein
Immunogen Species	Gerbilliscus gambianus (Gambian gerbil) (Gerbilliscus kempfi gambiana)
Purity	>85% (SDS-PAGE)
Sequence	TIWLDLNMFL SLGVDCWIDN TRVVYNRSSG RVSNA PGVEI RVPGF GKTY S VEYLDDNKLA EYMHTLVQNL VNNGYVRDET VRAAPYDWRL EPSQQDDYYQ KLAGLIEEMY AAYGKPVFLI GHSLGCLHVL YFLLRQGIPI MSSIKLREEQ RITTTSPWMF PDRDVPEDH VFISTPEFNY TGQDFERFFS DLHFEEGWYM WLQSRDLLAG LPAPGVDVYC LYGVGLPTPH TYIYDHNFPY KDPVAALYED GDDTVATRST ELCGQWQGRQ SQPVHLLPMN GTEHLNMVFS NKT
Source	Yeast
Target Names	LCAT
Protein Names	Recommended name: Phosphatidylcholine-sterol acyltransferase EC= 2.3.1.43 Alternative name(s): Lecithin-cholesterol acyltransferase Phospholipid- cholesterol acyltransferase
Expression Region	1-293
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.