



HLA-DPA1 Antibody, HRP conjugated

Product Code	CSB-PA14819B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20036
Immunogen	Recombinant Human HLA class II histocompatibility antigen, DP alpha 1 chain protein (29-222AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HLA class II histocompatibility antigen, DP alpha 1 chain (DP(W3)) (DP(W4)) (HLA-SB alpha chain) (MHC class II DP3-alpha) (MHC class II DPA1), HLA-DPA1, HLA-DP1A HLASB
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	HLA-DPA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.