





KRT18 Antibody, Biotin conjugated

Product Code	CSB-PA10629D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05783
Immunogen	Recombinant Human Keratin, type I cytoskeletal 18 protein (102-396AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin, type I cytoskeletal 18 (Cell proliferation-inducing gene 46 protein) (Cytokeratin-18) (CK-18) (Keratin-18) (K18), KRT18, CYK18
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	KRT18

Image



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA10629D0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA10629D0Rb at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.