



Actb Antibody, Biotin conjugated

Product Code	CSB-PA07675D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60710
Immunogen	Recombinant Mouse Actin, cytoplasmic 1 protein (11-375AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Actin, cytoplasmic 1 (Beta-actin) [Cleaved into: Actin, cytoplasmic 1, N-terminally processed], Actb
Immunogen Species	Mus musculus (Mouse)
Research Area	Isotype/Loading Controls
Target Names	Actb
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.