



MAOB Antibody, Biotin conjugated

Product Code	CSB-PA013399LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27338
Immunogen	Recombinant Human Amine oxidase [flavin-containing] B protein (2-489AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Amine oxidase [flavin-containing] B (EC 1.4.3.4) (Monoamine oxidase type B) (MAO-B), MAOB
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MAOB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.