



# Alpl Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA17039C0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P09242
<b>Immunogen</b>	Recombinant Mouse Alkaline phosphatase, tissue-nonspecific isozyme protein (20-503AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Alkaline phosphatase, tissue-nonspecific isozyme (AP-TNAP) (TNSALP) (EC 3.1.3.1) (Alkaline phosphatase 2) (Alkaline phosphatase liver/bone/kidney isozyme), Alpl, Akp-2 Akp2
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	Alpl