




CXCL13 Antibody

Product Code	CSB-PA06139A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43927
Immunogen	Recombinant Human C-X-C motif chemokine 13 protein (23-95AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC;Recommended dilution:IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.1mol/L NaCl, 15mmol/L NaN3, PH7.2
Purification Method	Caprylic Acid Ammonium Sulfate Precipitation purified
Isotype	IgG
Clonality	Polyclonal
Alias	Angie,B cell-attracting chemokine 1,B lymphocyte chemoattractant,CXC chemokine BLC,Small-inducible cytokine B13,CXCL13,BCA1, BLC, SCYB13
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CXCL13

Image



Immunohistochemistry of paraffin-embedded human liver using CSB-PA06139A0Rb at dilution 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.