



PAX2 Antibody

Product Code	CSB-PA11564A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02962
Immunogen	Recombinant Human Paired box protein Pax-2 protein (1-363AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Paired box protein Pax-2, PAX2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PAX2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.