



fliC Antibody, Biotin conjugated

Product Code	CSB-PA09734D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04949
Immunogen	Recombinant Escherichia coli Flagellin protein (2-498AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Flagellin, fliC, flaF hag
Immunogen Species	Escherichia coli (strain K12)
Research Area	Immunology
Target Names	fliC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.