



AQP1 Antibody, Biotin conjugated

Product Code	CSB-PA16085D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29972
Immunogen	Recombinant Human Aquaporin-1 protein (220-269AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aquaporin-1 (AQP-1) (Aquaporin-CHIP) (Urine water channel) (Water channel protein for red blood cells and kidney proximal tubule), AQP1, CHIP28
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AQP1