



AQP1 Antibody, HRP conjugated

Product Code	CSB-PA16085B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29972
Immunogen	Recombinant Human Aquaporin-1 protein (220-269AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aquaporin-1 (AQP-1) (Aquaporin-CHIP) (Urine water channel) (Water channel protein for red blood cells and kidney proximal tubule), AQP1, CHIP28
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AQP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.