



viaA Antibody, Biotin conjugated

Product Code	CSB-PA16519D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8Z2R0
Immunogen	Recombinant Salmonella typhi Protein ViaA protein (3-476AA)
Raised In	Rabbit
Species Reactivity	Salmonella typhi
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein ViaA (VWA domain protein interacting with AAA ATPase), viaA
Immunogen Species	Salmonella typhi
Research Area	Others
Target Names	viaA