



viaA Antibody, HRP conjugated

Product Code	CSB-PA16519B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8Z2R0
Immunogen	Recombinant Salmonella typhi Protein ViaA protein (3-476AA)
Raised In	Rabbit
Species Reactivity	Salmonella typhi
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein ViaA (VWA domain protein interacting with AAA ATPase), viaA
Immunogen Species	Salmonella typhi
Research Area	Others
Target Names	viaA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.