



# Recombinant Differentially expressed in FDCP 8 homolog (F20H11.1)

<b>Product Code</b>	CSB-BP006651CXY
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O01738
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSSWCSSEDA HMQSSTPSTR SRKSSWLISK IDENEIDKEV QEMMDDWRRS NEAENHVECI KAALEYCRDK LKELLIDEED TEMDGKETKT KPINIVSNTE FPEDRDVLGN LVNCLDRQQK AARQQMYFDK IVQLGLERQA LQDEAVVSTS SPLDECRAEG HEFVMQPVRG GHNPCCEVCM HTIWRLVQWW RRCRVCGMRA HDKCAEEVKR VCAGVLSTRA KFELNTNLCE ERSLAEQEYQ CAECTAPICF DGVAEQEARL CDYSGELFCP NCHWNDVWSI PARIVHNLDS TPRPVCRAVK QLLAIDHRP LIDINESTLS LIKFHKELRR VNELRRNLL MKCYFVSCRT ARRLRILQYL NAHSHFVDNS VMYSPKELRE LCDGTLPLDL EQIHTVFRKH IEEECETCAG NGFFCELCD INVDQKNKIL YPFTENTRSC ATCLAVYHKK CFERKSLNCP RCERRRKRTE IPKTLSSSDE KEAVRN
<b>Source</b>	Baculovirus
<b>Target Names</b>	F20H11.1
<b>Protein Names</b>	Recommended name: Differentially expressed in FDCP 8 homolog
<b>Expression Region</b>	1-486
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.