



Recombinant Differentially expressed in FDCP 8 homolog (F20H11.1)

Product Code	CSB-EP006651CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O01738
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MSSWCSSEDA HMQSSTPSTR SRKSSWLISK IDENEIDKEV QEMMDDWRRS NEAENHVECI KAALEYCRDK LKELLIDEED TEMDGKETKT KPINIVSNTE FPEDRDVLGN LVNCLDRQQK AARQQMYFDK IVQLGLERQA LQDEAVVSTS SPLDECRAEG HEFVMQPVRG GHNPCCEVCM HTIWRLVQWW RRCRVCGMRA HDKCAEEVKR VCAGVLSTRA KFELNTNLCE ERSLAEQEYQ CAECTAPICF DGVAEQEARL CDYSGELFCP NCHWNDVWSI PARIVHNLDS TPRPVCRAVK QLLAIDHRP LIDINESTLS LIKFHKELRR VNELRRNLL MKCYFVSCRT ARRLRILQYL NAHSHFVDNS VMYSPKELRE LCDGTLPLDL EQIHTVFRKH IEEECETCAG NGFFCELDD INVQKKNIL YPFTENTRSC ATCLAVYHKK CFERKSLNCP RCERRRKRTE IPKTLSSSDE KEAVRN
Source	E.coli
Target Names	F20H11.1
Protein Names	Recommended name: Differentially expressed in FDCP 8 homolog
Expression Region	1-486
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.