



VEGFD Antibody, FITC conjugated

Product Code	CSB-PA07554C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43915
Immunogen	Recombinant Human Vascular endothelial growth factor D protein (89-205AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vascular endothelial growth factor D (VEGF-D) (c-Fos-induced growth factor) (FIGF), VEGFD, FIGF
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	VEGFD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.