




IL-33 Antibody

Product Code	CSB-PA011656DA01PI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	M5B263
Immunogen	Recombinant Pig Interleukin-33 protein (123-276AA)
Raised In	Rabbit
Species Reactivity	Pig, Mouse, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	IL-33
Immunogen Species	Pig
Research Area	Others
Target Names	IL-33

Image



Western blot
 All lanes: IL-33 antibody at 2µg/ml
 Lane 1: HeLa whole cell lysate
 Lane 2: Mouse lung tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 32 kDa
 Observed band size: 32 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.