





Avpr1a Antibody

Product Code	CSB-PA15624A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30560
Immunogen	Recombinant Rat Vasopressin V1a receptor protein (7-52AA)
Raised In	Rabbit
Species Reactivity	Rat, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Neuroscience
Target Names	Avpr1a

Image



Western blot
 All lanes: Avpr1a antibody at 2µg/ml + Mouse spleen tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 48 kDa
 Observed band size: 48 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA15624A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.