



H2AFJ Antibody, FITC conjugated

Product Code	CSB-PA010094DC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BTM1
Immunogen	Recombinant Human Histone H2A.J protein (2-129AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Histone H2A.J (H2a/j), H2AFJ
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	H2AFJ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.