



AKT1 Antibody, Biotin conjugated

Product Code	CSB-PA15905D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31749
Immunogen	Recombinant Human RAC-alpha serine/threonine-protein kinase protein (1-480AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RAC-alpha serine/threonine-protein kinase (EC 2.7.11.1) (Protein kinase B) (PKB) (Protein kinase B alpha) (PKB alpha) (Proto-oncogene c-Akt) (RAC-PK-alpha), AKT1, PKB RAC
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	AKT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.