



Thbs1 Antibody, FITC conjugated

Product Code	CSB-PA023487LC01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35441
Immunogen	Recombinant Mouse Thrombospondin-1 protein (19-350AA)
Raised In	Rabbit
Species Reactivity	Mouse
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thrombospondin-1, Thbs1, Tsp1
Immunogen Species	Mus musculus (Mouse)
Research Area	Neuroscience
Target Names	Thbs1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.