



PIGR Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA017981LC01BO |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P81265 |
| Immunogen | Recombinant Bovine Polymeric immunoglobulin receptor protein (399-599AA) |
| Raised In | Rabbit |
| Species Reactivity | Bovine |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Polymeric immunoglobulin receptor (PIgR) (Poly-Ig receptor) [Cleaved into: Secretory component], PIGR |
| Immunogen Species | Bos taurus (Bovine) |
| Research Area | Others |
| Target Names | PIGR |