




PIGR Antibody

Product Code	CSB-PA017981LA01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P81265
Immunogen	Recombinant Bovine Polymeric immunoglobulin receptor protein (399-599AA)
Raised In	Rabbit
Species Reactivity	Bovine
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polymeric immunoglobulin receptor (PIgR) (Poly-Ig receptor) [Cleaved into: Secretory component], PIGR
Immunogen Species	Bos taurus (Bovine)
Research Area	Others
Target Names	PIGR

Image



Western blot

All lanes: Polymeric immunoglobulin receptor antibody at 2µg/ml + HL60 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 83, 59 kDa

Observed band size: 83 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.