



CST3 Antibody, HRP conjugated

Product Code	CSB-PA06197B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01034
Immunogen	Recombinant Human Cystatin-C protein (27-146AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cystatin-C (Cystatin-3) (Gamma-trace) (Neuroendocrine basic polypeptide) (Post-gamma-globulin), CST3
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	CST3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.