



# Efna1 Antibody

Product Code	CSB-PA007460LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P97553
Immunogen	Recombinant Rat Ephrin-A1 protein (18-182AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ephrin-A1 (EPH-related receptor tyrosine kinase ligand 1) (LERK-1) (Immediate early response protein B61) [Cleaved into: Ephrin-A1, secreted form], Efna1, Epgl1 Lerk1
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cardiovascular
Target Names	Efna1