



# RBP1 Antibody

<b>Product Code</b>	CSB-PA019480GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P09455
<b>Immunogen</b>	Human RBP1 fusion protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA,WB,IHC
<b>Storage Buffer</b>	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen affinity purification
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	CRABP-I; CRBP; CRBP1; CRBPI; RBPC
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	RBP1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.