



Ccl9 Antibody

Product Code	CSB-PA09119A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51670
Immunogen	Recombinant Mouse C-C motif chemokine 9 protein (22-114AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.1mol/L NaCl, 15mmol/L NaN ₃ , PH7.2
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-C motif chemokine 9 (CCF18) (Macrophage inflammatory protein 1-gamma) (MIP-1-gamma) (Macrophage inflammatory protein-related protein 2) (MRP-2) (Small-inducible cytokine A9) [Cleaved into: CCL9(29-101); CCL9(30-101); CCL9(31-101)], Ccl9, Mrp2 Scya10 Scya9
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Ccl9