



Ccl9 Antibody

| | |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-PA09119A0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P51670 |
| Immunogen | Recombinant Mouse C-C motif chemokine 9 protein (22-114AA) |
| Raised In | Rabbit |
| Species Reactivity | Mouse |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.1mol/L NaCl, 15mmol/L NaN ₃ , PH7.2 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | C-C motif chemokine 9 (CCF18) (Macrophage inflammatory protein 1-gamma) (MIP-1-gamma) (Macrophage inflammatory protein-related protein 2) (MRP-2) (Small-inducible cytokine A9) [Cleaved into: CCL9(29-101); CCL9(30-101); CCL9(31-101)], Ccl9, Mrp2 Scya10 Scya9 |
| Immunogen Species | Mus musculus (Mouse) |
| Research Area | Immunology |
| Target Names | Ccl9 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |