



metC Antibody

Product Code	CSB-PA14479A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06721
Immunogen	Recombinant Escherichia coli Cystathionine beta-lyase MetC protein (1-392AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cystathionine beta-lyase MetC (CBL) (CL) (EC 4.4.1.8) (Beta-cystathionase) (Cysteine desulfhydrase) (CD) (EC 4.4.1.28) (Cysteine lyase), metC
Immunogen Species	Escherichia coli
Research Area	Signal Transduction
Target Names	metC