




# Avpr1b Antibody

<b>Product Code</b>	CSB-PA15639A0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P48974
<b>Immunogen</b>	Recombinant Rat Vasopressin V1b receptor protein (343-425AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Rat
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Vasopressin V1b receptor (V1bR) (AVPR V1b) (AVPR V3) (Antidiuretic hormone receptor 1b) (Vasopressin V3 receptor), Avpr1b
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Research Area</b>	Neuroscience
<b>Target Names</b>	Avpr1b

## Image



### Western blot

All lanes: Vasopressin V1b receptor antibody at 2µg/ml

Lane 1: EC109 whole cell lysate

Lane 2: 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 47 kDa

Observed band size: 47 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.