




Ctsk Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA006192LA01RA |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O35186 |
| Immunogen | Recombinant Rat Cathepsin K protein (115-329AA) |
| Raised In | Rabbit |
| Species Reactivity | Rat |
| Tested Applications | ELISA, IHC; Recommended dilution: IHC:1:10-1:100 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Cathepsin K (EC 3.4.22.38), Ctsk |
| Immunogen Species | Rattus norvegicus (Rat) |
| Research Area | Cardiovascular |
| Target Names | Ctsk |

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA006192LA01RA at dilution of 1:20

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.