




Ctsk Antibody

Product Code	CSB-PA006192LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O35186
Immunogen	Recombinant Rat Cathepsin K protein (115-329AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:10-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cathepsin K (EC 3.4.22.38), Ctsk
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cardiovascular
Target Names	Ctsk

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA006192LA01RA at dilution of 1:20