



gadA Antibody, Biotin conjugated

Product Code	CSB-PA302526HD01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P69908
Immunogen	Recombinant Escherichia coli Glutamate decarboxylase alpha protein (1-466AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamate decarboxylase alpha (GAD-alpha) (EC 4.1.1.15), gadA, gadS
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	gadA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.